

CHARACTERIZATION OF QUANTUM DOTS CONJUGATES WITH ANTIBODIES BY CAPILLARY ELECTROPHORESIS




Ivona Voracova, Marcela Liskova, Karel Kleparnik, Vera Hezinova, Jan Prikryl, Frantisek Foret
 Institute of Analytical Chemistry of the ASCR, v.v.i., Veveří 97, 602 00 Brno, Czech Republic, voracova@iach.cz



Preparation



1st step: preparation of hydrogen telluride



2nd step: quantum dots formation

$\text{CdCl}_2 + \text{NaHTe} + \text{MPA}$ or $\text{MA} + \text{heating}$

MPA: 3-mercaptopropionic acid: $\text{HS-CH}_2\text{-CH}_2\text{-COOH}$

MA: 2-mercaptoethylamin: $\text{HS-CH}_2\text{-CH}_2\text{-NH}_2$


3rd step: coating

CdS: $\text{CdTe} + \text{CdCl}_2 + \text{Na}_2\text{S} + \text{MPA}$ or $\text{MA} + \text{heating}$

ZnS: $\text{CdTe/CdS} + \text{ZnCl}_2 + \text{Na}_2\text{S} + \text{MPA}$ or $\text{MA} + \text{heating}$




Luminescence spectra



- wide excitation spectra with maximum at 469 nm

- narrow emission spectra with maximum at 600 nm
 - bandwidth 58 nm at half height

Emission spectra



Quantum yield


Coating	Ligand	
	MPA	MA
-	9.66	21.08
CdS	18.18	25.54
CdS/ZnS	9.08	27.96

$$QY_{QD} = QY_S \frac{A_S I_{QD}}{I_S A_{QD}}$$


S – standard
 (fluorescein QY = 0.90*)

* Demas J.N, Crosby G.A., J Phys. Chem., 75, 991 (1971)

Luminescence lifetimes




Electrophoresis in replaceable sieving media




PVA coated capillary 20/30 cm, i.d. 75µm, separation buffer 3% LPA 10 MDa in 50 mM TRIS/TAPS buffer, pH = 9, QD 2.8 and 3.7 nm (525 and 610 nm 1:1), injection time 10 s, separation voltage 3 kV

Conjugation using zero-cross linkers



- formation of peptide bond
- carboxylic group of MPA on the QD surface
- amino group from antibody
- catalyzers: EDC (carbodiimide)
HNS (succinimide)


Capillary electrophoresis of conjugates



Separation conditions: fused silica capillary 15/25cm, separation voltage 6 kV, injection time 15 s, separation buffer 100 mM TRIS/TAPS, pH = 8.3


Capillary electrophoresis of conjugates

Apoptotic cells of mouse duodenum labeled by QD-annexin conjugate




Separation conditions: fused silica capillary 15/25cm, separation voltage 6 kV, injection time 15 s, separation buffer 100 mM TRIS/TAPS, pH = 8.3

Capillary electrophoresis of conjugates




Separation conditions: fused silica capillary 15/25 cm, separation voltage 6 kV, injection time 15s, separation buffer 100 mM TRIS/TAPS

Capillary electrophoresis of conjugates





Separation conditions: fused silica capillary 15/25 cm, separation voltage 6 kV, injection time 6s, separation buffer 100 mM TRIS/TAPS

Conjugation of magnetic beads with Anti-ovalbumin



Capillary electrophoresis of conjugates

- conjugated with antibody against CD3 protein – a membrane protein specific for T-lymphocytes while B-lymphocytes do not contain it



Separation conditions: fused silica capillary 12/20 cm, separation voltage 6 kV, injection time 15 s, separation buffer 50 mM CAPS, pH = 11.2

Conclusions

CdTe quantum dots

- fluorescence lifetimes depend on particle size
- small selectivity of size separation using sieving media
- conjugation with antibody and macrocyclic ligand
- capillary electrophoresis luminescence immunoassay
- selective cell labeling

Acknowledgment

HRTEM

Dr. Mariana Klementová – Institute of Inorganic Chemistry AS CR

Fluorescence and fluorescence lifetimes measurement

Mgr. Petr Táborský, Ph.D. and doc. Jan Preisler, Ph.D.
 – Department of Chemistry, Faculty of Science, Masaryk University

Financial support

GA AV - KAN400310651 and KJB400310709
 GACR - GA203/08/1680
 MSMT - LC06023
 AV0Z40310501